

# clinell®

## Spill Wipes Kit

---

Complete kit for clearing body  
fluid spills, from the UK's no.1  
hospital disinfectant brand



# Clinell Spill Wipes Kit

Absorbs an entire body fluid spill for safe removal in seconds.

**Kills microorganisms associated with body fluid spills**

Proven efficacy against Hepatitis B & C, *Clostridioides difficile*, norovirus and more, killing all known microorganisms including non-enveloped viruses and bacterial endospores. Conforms to: EN17126, EN13727, EN14348, EN14476 & EN13624.

**Quicker alternative to traditional spill kits**

Replaces traditional chlorine-based spill kits, which are time consuming and hazardous to use. Clinell Spill Wipes are simple to use, helping to increase compliance and reduce the burden of pathogen transmission.

**Based on patented technology**

Harnessing the formulation technology of Clinell Peracetic Acid Wipes\*, the pad generates powerful oxidising agents including peracetic acid, which is safe for use on blood, faeces and urine - unlike chlorine-based disinfectants, which form toxic gas on contact with urine.

**Soak spill, wipe surface, dispose safely**

As a complete kit, Clinell Spill Wipes provide an efficient solution for spills up to 350ml from start to finish. After use, the soiled pad, wipes and PPE worn by the user can be placed into the re-sealable bag and safely disposed of with the appropriate waste channel.

## What’s in a Clinell Spill Wipes Kit?



One super absorbent, peracetic acid generating pad



Two large detergent-disinfectant wipes



Re-sealable bag for safe and easy waste disposal

# Body fluid spills

Designed to make dealing with blood, urine and faecal spills quicker and easier.

Traditionally, body fluid spills are decontaminated with chlorine before safe removal. This lengthy process can take at least 10 minutes. One spill typically requires multiple steps and products, making compliance a challenge to everyday users.

**Removing a spill with chlorine-based spill kits**

- Place barrier agent around the edges of the spill
  - Pour 10,000ppm chlorine granules over the spill
  - Allow time to neutralise pathogens
- Remove granules with scoop or dustpan
  - Dispose of in biohazard bag
  - Wipe the area with absorbent towels to remove remaining residue

Chlorine-based spill kits can expose users to hazardous chemicals and can’t be used on urine – as they generate carcinogenic toxic fumes on contact. This means that users need two different kits: one for blood and faeces, and one for urine. With Clinell Spill Wipes, one kit can be used for all body fluid spills.



Product	Unit of issue	Order code	NHSCC
Spill Wipes	Single unit	CSW1	VJT268
Wall Mounted Dispenser	(for CSW1)	CSW1D	-





**GAMA Healthcare Ltd.,**

The Maylands Building, Maylands Avenue,  
Hemel Hempstead, HP2 7TG, UK.

T: +44 (0)20 7993 0030

E: [info@gamahealthcare.com](mailto:info@gamahealthcare.com)

[www.gamahealthcare.com](http://www.gamahealthcare.com)

**Caution:**

Use biocides safely. Always read the label  
and product information before use.

JBN220409